



volume	10uL, 20uL, 100uL, 200uL, 300uL, 1000uL, 1250uL
Tip Color	Clear / Selective
Filter	Filtered / Non Filtered
Sterilization	Sterilize / Selective
Material	Polypropylene
Rack Color	Different colors
Test Passed	Free of DNase, RNase, and Human DAN
Certificate	ISO13485/ CE
Compatible Pipettor	Eppendorf Pipettor Gilson Pipettor Thermo Pipettor Brand Pipettor YongYue Pipettor

比例 SCALE	1:1
单位 UNIT	MM
SHEET	OF

 永岳医疗科技(昆山)有限公司 Yong Yue Medical Technology (Kunshan) Co., Ltd.	
图号 DRAWING NO	
图名 DRAWING NAME	Filter tips 10uL, 20uL, 100uL, 200uL, 300uL, 1000uL, 1250uL
设计 DRAWING BY	Cheng
审核 DESIGNED BY	Yao
核准 APPROVED BY	Yao